Active.valve TF

6P20/W




Instructions for use

MD

Active.valve TF

DE	Gebrauchsanweisung Oberschenkel Ausstoßventil, wasserfest.....	5
EN	Instructions for use AK exhaust valve, waterproof.....	7
ES	Instrucciones de uso Válvula de descarga para prótesis transfemoral, resistente al agua....	9
FR	Notice d'utilisation Valve d'expulsion fémorale, étanche à l'eau.....	12
IT	Istruzioni per l'uso Valvola di espulsione femorale, impermeabile	14
PL	Instrukcja używania Zawór wypychający uda, wodoodporny	16
CS	Návod k použití Vytěšňovací ventil na stehenní protézy, vodotěsný	19

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■ Nutzung

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EN Instructions for use

■ Use

The **6P20/W Active.valve TF** (hereinafter referred to as "the product") is approved for **activity levels 1 to 4** and may be used only for light sporting activities (e.g. light jogging < 7 km/h). The product must not be exposed to overloading, i.e. impacts, shocks or vibrations.

The indication-appropriate use must be assessed by a prosthetic specialist.

The product may be fitted only by a prosthetic specialist.

This product is intended for single patient use only.

Please note that this product is tested and approved only in combination with other Streifeneder modular components.

■ Useful life

The maximum useful life of the product is limited to 36 months.

■ Intended Use

The product is exclusively intended for lower extremity prosthetic treatment.

■ Assembling

In the distal area, the product can be attached medially or laterally. With the enclosed placeholders (dummies) (**Fig. 1, A (deep-drawing dummy) and B (laminating dummy)**), which are used for deep-drawing or laminating, precise fitting in the different socket materials is possible.

An airtight socket is the basic requirement for use of the product. It is therefore necessary to manufacture the socket with the greatest care according to one of the methods specified in these instructions for use. We do not assume any liability for other socket constructions.

If socket adjustments or grinding tasks should be necessary, please seal the openings of the product on the inside and outside of the socket with polyethylene adhesive tape.

First, the positive cast is modelled in such a way that the deep-drawing dummy has a level surface. Ensure that there is no gap between the deep-drawing dummy and the positive cast.

Pull a vacuum stockinette/fine stockinette over the positive cast.

The deep-drawing dummy is fixed to the positive cast using the screw supplied (Fig. 2).

■ Procedure in the deep-drawing process

Deep-draw the test socket or the elastic inner socket (Fig. 3). Then grind the deep-drawing dummy flush and unscrew the screw (Fig. 4).

For plastic inner sockets that are laminated over, the lamination dummy is fixed on the flat ground support surface of the deep-drawing dummy with the screw (Figs. 5, 6 and 7).

■ Procedure in the lamination process

For designs with an inner shaft, see above for use of the lamination dummy.

For manufacturing a resin socket without inner socket, only attach the deep-drawing dummy to the positive cast (Fig. 2).

- Pull the soaked PVA film tube over the dummy/dummies and drain
- Apply reinforcement materials depending on the expected load
- Pull over the soaked PVA film tube and carry out the casting process

Grind the dummy flat after lamination, unscrew the screw and separate the positive cast from the socket (Fig. 5).

The surfaces must be ground flat and parallel to ensure tightness.

■ Replacing the spring

To remove the spring, open the cover by turning it anti-clockwise using an Allen key (size 4 mm) (Fig. 8).

The spring (from 6P20/WE99) must be aligned centrically in the recess of the cover (Fig. 9). Screw the cover back into the valve body using the Allen key while making sure that the spring is aligned.

■ Replacing the O-ring

If the product does not sufficiently maintain the vacuum, a worn O-ring may be the cause.

To replace the old O-ring, use a small slotted screwdriver for example to remove it. The new O-ring (from 6P20/WE99) is then carefully placed in the intended groove (Fig. 10, A).

Alternatively, the complete valve body (6P20/WE98, Fig. 10) can be replaced.

■ Technical Notes

It is essential to ensure that the product is not exposed to chemicals and substances such as acids, dust or sand!

The corrosion resistance of this product is limited and depends on the component with the lowest corrosion resistance used in the prosthesis structure, which determines the maximum corrosion resistance of the entire structure.

The prosthesis must be rinsed thoroughly with fresh water every time it comes into contact with water. After rinsing, the prosthesis should be dried.

Failure to observe these instructions will make any claims against Streifeneder ortho.production GmbH null and void.

Please inform your patient accordingly!

■ Function and application

The valve body can be removed from the prosthetic socket by pressing the locking slides together. For example, for threading through a donning aid.

To release the vacuum in the prosthetic socket, press the aluminium button in the centre.

■ Maintenance and Service

Generally, all movable parts must be inspected visually and for functionality every 6 months; all non-movable structural parts every 12 months.

■ Declaration of Conformity

Streifeneder ortho.production GmbH declares as manufacturer with sole responsibility that the product complies with Regulation (EU) 2017/745 for medical devices (MDR).

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© copyright Streifeneder ortho.production