Shuttle-Lock

15S2/W





Instructions for use

MD

Shuttle-Lock

DE	Gebrauchsanweisung	4
EN	Instructions for use	6
ES	Instrucciones de uso	8
FR	Notice d'utilisation	11
IT	Istruzioni per l'uso	13
PL	Instrukcja używania	16
CS	Návod k použití	18

■ Nutzung

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EN Instructions for use

■ Use

The **15S2/W Shuttle-Lock** (hereinafter referred to as “the product”) is approved for **activity levels 1 to 4** and may be used only for light sporting activities (e.g. light jogging < 7 km/h). The product must not be exposed to overloading, i.e. impacts, shocks or vibrations.

The indication-appropriate use must be assessed by a prosthetic specialist.

The product may be fitted only by a prosthetic specialist.

This product is intended for single patient use only.

Please note that this product is tested and approved only in combination with other Streifeneder modular components.

■ Useful life

The maximum useful life of the product is limited to 36 months.

■ Intended use

The product is exclusively intended for lower extremity prosthetic treatment and can be used in fresh, salt and chlorinated water.

■ Assembling

The specified build height of 23 mm (housing, **A**) or 56 mm (housing and lamination protection cone, **B**) enables the other modular parts to be determined when setting up the prosthesis (**Fig. 1**).

If necessary, unscrew the locking unit anticlockwise from the Shuttle-Lock (**Fig. 2**).

Reduce the distal end of the stump on the positive cast to a surface with a diameter of 47 mm, perpendicular to the direction of the stump (Fig. 3).

■ Reinforcing the product

The reinforcement instructions described here are the minimum standard. Additional layers must be added, depending on the expected loads.

Fix the 10S1/E6L lamination protection cone to the positive cast using the enclosed nails.

Carefully grease all threads and/or surfaces that should not be penetrated by resin, or insulate with 113P7 plasticine.

- Cover with the soaked 119P3 PVA lamination bag.
- Screw in the short locking unit dummy (white) and mount the lock with all other dummies on the 10S1/E6L.
- Tighten the dummy hand-tight so that the O-ring (A) is in contact but not squeezed (Fig. 4).
- Cover with 2 layers of 95P17 polyglass stockinette.
- Apply 1 layer of 95P8 carbon fibre tape in the ML direction.
- Apply 1 layer of 95P8 carbon fibre tape in the AP direction.
- Apply 1 layer of 95P8 carbon fibre tape in the notches.
- Cover with 2 layers of 95P9 carbon fibre braided sleeve and tie in the notches provided on the product.
- Apply further layers, depending on the expected loads.
- Cover with 2 layers of 95P17 polyglass stockinette.
- Cover with the soaked 119P3 PVA lamination bag and tie proximally. Carry out the casting process with 112P120 AcryliPure lamination resin or 112P122 AcryliPure lamination resin C.

After curing, remove the 119P3 PVA lamination bag and roughen the socket.

For foaming, use 112P25 rigid foam H700 and the long locking unit dummy (white).

For further reinforcement, see the instructions for use of the components used.

Do not remove the locking unit dummy until all sanding work has been completed.

Subsequent drilling through the pin guide is not permitted!

The release button of the locking unit must always be accessible from the outside.

■ Technical notes

The locking unit can be opened by turning anticlockwise (right-handed thread).

Tighten the locking unit to **5 Nm** using the calibrated 173P10 torque spanner.

A socket wrench (size 16) and, if applicable, an adapter for torque spanners from 1/4" to 3/8" is/are required.

It is essential to ensure that the product is not exposed to corrosive chemicals and substances such as acids, dust or sand!

The corrosion resistance of this product is limited and depends on the component with the lowest corrosion resistance used in the prosthesis structure, which determines the maximum corrosion resistance of the entire structure.

The prosthesis must be rinsed thoroughly with fresh water every time it comes into contact with water. After rinsing, the prosthesis should be dried.

Failure to observe these instructions will make any claims against Streifeneder ortho. production GmbH null and void.

Please inform your patient accordingly!

■ Replacing the push button

When replacing the push button, first clean the thread and then secure the new push button with superglue.

■ Maintenance and Service

Generally, all movable parts must be inspected visually and for functionality every 6 months; all non-movable structural parts every 12 months. If the product is also used for light sporting activities, the inspection intervals must be reduced to a maximum of 3 months. If the product is also used in water, the maintenance intervals must be reduced to a maximum of 3 months.

In the event of malfunctions, the product must be sent in.

Please note that repairs to the product may be carried out only by Streifeneder ortho. production GmbH.

■ Declaration of conformity

Streifeneder ortho.production GmbH declares as manufacturer with sole responsibility that the product complies with Regulation (EU) 2017/745 for medical devices (MDR).

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© copyright Streifeneder ortho.production
Upozorňujeme, že opravy výrobku smí provádět pouze společnost Streifeneder ortho.production GmbH.