KINEGEN.flow

3A1900

Instructions for use


MD

KINEGEN.flow


DE	Gebrauchsanweisung Polyzentrisches Kniegelenk	7
EN	Instructions for use Polycentric knee joint.....	11
ES	Instrucciones de uso Rodilla policéntrica.....	16
FR	Notice d'utilisation Articulation de genou polycentrique.....	20
EL	Οδηγίες χρήσης Πολυκεντρική άρθρωση γόνατος.....	25
IT	Istruzioni per l'uso Articolazione policentrica del ginocchio	30
PL	Instrukcja używania Policentryczny przegub kolanowy	35
CS	Návod k použití Polycentrický kolenní kloub.....	39
SV	Bruksanvisning Polycentrisk knäled	43
DA	Brugsanvisning Polycentrisk knæled	48

PT	Instruções de utilização Articulação policêntrica do joelho	52
CP	Упутство за употребу Полицентрични зглоб колена.....	57
BS	Uputstvo za upotrebu Policentrični zglob koljena	61
RO	Instrucțiuni de utilizare Articulație de genunchi policentrică.....	66
SK	Návod na použitie Polycentrický kolenný kĺb	70
BG	Инструкции за употреба Полицентрична колянна става	75
FI	Käyttöohje Monikeskinen polvinivel.....	80
LT	Naudojimo instrukcija Daugiacentris kelio sąnarys	84
ET	Kasutusjuhend Polütsentriline põlveliiges.....	88


1


2


3


4




5

6

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■ **Nutzung**

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EN Instructions for use

■ Use

The **3A1900 KINEGEN.flow** (hereinafter referred to as “the product”) is approved for **activity levels 2 to 4** and for a maximum body weight, including carried and supported loads, of up to **125kg**, and may be used only for light sporting activities (e.g. light jogging < 7 km/h). The product must not be exposed to overload, i.e. impacts, shocks or vibrations.

The approved maximum weight load of this product is limited and depends on the component with the lowest weight limit used in the prosthesis structure, which determines the maximum weight load of the entire structure.

The indication-appropriate use must be assessed by a prosthetic specialist.
The product may be fitted only by a prosthetic specialist.
This product is intended for single patient use only.

Please note that this product is tested and approved only in combination with other Streifeneder and College Park modular components.

■ Useful life

The maximum useful life of the product is limited to 36 months.

■ Intended use

The product is exclusively intended for lower extremity prosthetic treatment and can be used in fresh, salt and chlorinated water.

■ Assembling

The arrangement of the knee joint axes ensures that, in the extended position, the position of pivot point (D) is above the knee joint and behind the load line (Fig. 1).

The stance phase stability is thus determined exclusively via the prosthetic alignment.

The product's properties can only be optimally utilised with a correctly aligned prosthesis. To ensure correct alignment, the position of the residual limb must be taken into account when positioning the socket connection.

Do not touch the joint mechanism when flexing the knee joint – **risk of injury from pinching fingers!**

■ Basic set-up

The prosthesis should be set up so that in the frontal view (Fig. 2) the plumb line (a) runs downwards from the centre of the socket over the centre of the proximal connection and through the distal adjustment core.

In the sagittal view (Fig. 3), the set-up reference line (b) runs from the greater trochanter (T) across the anterior upper joint axis to the transition between the posterior third and central third of the prosthetic foot, and the heel height is set to approx. 5 mm. The set-up reference point is indicated in the area of the anterior upper joint axis in the cover (Fig. 4, A).

The individual circumstances of the patient (e.g., the flexion contracture) must be taken into account in the further course.

Fine adjustments are made during the dynamic fitting.

If residual limb flexion is not taken into account, the centre of rotation of the joint may be too far forward, which can lead to functional disorders and even loss of stability.

■ Flexion limitation

The product enables knee flexion of 155°. However, this flexion angle can only be achieved with a correctly set up socket/prosthesis structure. The clamps of lamination anchors must be located medially or laterally (3-arm anchor), or diagonally anterior (4-arm anchor) on the socket. If the clamp is located posteriorly, the flexion angle of the joint may not be fully utilised or the joint may be damaged.

The product has an integrated flexion limitation. It is essential to ensure that the socket is not used as a limitation. This can lead to defects and/or complete failure of the prosthesis.

■ Technical notes

It is essential to ensure that this product is not exposed to corrosive chemicals and substances such as acids, dust or sand!

The corrosion resistance of this product is limited and depends on the component with the lowest corrosion resistance used in the prosthesis structure, which determines the maximum corrosion resistance of the entire structure.

The prosthesis must be rinsed thoroughly with fresh water every time it comes into contact with water. After rinsing, the prosthesis should be dried.

Failure to observe these instructions will make any claims against Streifeneder ortho.production GmbH null and void.

Please inform your patient accordingly!

To eliminate sounds in the cosmetic soft foam cover, use silicone release spray, not talcum powder. Using talcum powder increases the risk of the patient falling, as talcum powder withdraws grease from the mechanical components, possibly leading to malfunctions. If the product is used with talcum powder, all claims against Streifeneder ortho.production GmbH are null and void. Spray the silicone spray directly onto the friction surfaces of the cosmetic soft foam cover.

Furthermore, when using a cosmetic soft foam cover, make sure that it does not get stuck between the moving parts of the knee joint, as this can lead to a fall.

■ Adjusting the stance-phase damping

The stance-phase damping mechanism is configured so that the joint chain stabilises as it compresses. Nevertheless, some patients may require an adjustment period before they become fully confident with the stance-phase damping.

To increase the pre-tensioning on the bumper and thus reduce the compression, loosen the clamping screw (**Fig. 5, K**) using an Allen key (size 4 mm), then turn the adjusting screw (**Fig. 6, J**) counterclockwise. An Allen key (size 2 mm) can be used for this for example. This results in a firmer heel response. As bumpers come in different hardness levels, the full adjustment range cannot be utilised with all bumpers. Excessive force may damage the tool or joint, even causing it to break.

For a softer heel strike, turn the adjusting screw (**Fig. 6, J**) clockwise. Ensure that the bumper is always sufficiently pre-tensioned. If there is play between the bumper and the adjustment nut, the adjusting screw should be turned approximately half a turn (three adjustment steps) counterclockwise after it makes contact with the bumper.

The bumper must always be pre-tensioned to ensure safe joint operation.

After the adjustment is completed, the locking screw (**Fig. 5, K**) must be tightened again to prevent unintentional adjustment of the damping.

For approximate adjustment of the damping behaviour, the bumper may be exchanged. In addition to the factory-fitted bumper with medium hardness (**M**), a softer bumper (**S**) and a harder bumper (**H**) are also available.

Depending on the patient's requirements, the table below can be used as a guide for selecting bumpers:

Hardness	Marking	Order number
soft	S	3A1900/E74
medium	M	3A1900/E75
hard	H	3A1900/E76

To replace the bumper, use an Allen key (size 4 mm) to loosen and remove the clamping screw (**Fig. 5, K**), turn the adjusting screw (**Fig. 6, J**) all the way down, and use a slotted screwdriver (size 1.8 mm) to completely unscrew the two locking screws (**Fig. 5, S**) on both sides of the bearing. The bearing can now be swung out to the rear (**Fig. 7**) and the bumper can be replaced. Then reassemble in reverse order.

Tighten the clamping screw (**Fig. 5, K**) to 2 Nm using the calibrated 173P11 torque spanner.

Tighten the locking screws (**Fig. 5, S**) hand-tight.

■ Adjusting the hydraulic swing phase control

The sequence for adjusting the swing phase control is as follows:

1. Using an Allen key (size 2 mm), reduce the flexion resistance by turning the valve screw (**Fig. 8, F**) anti-clockwise until the heel rise on the prosthesis side corresponds to the contralateral side. For less heel rise, turn the valve screw clockwise. Clicks and the scale on the valve screw enable precise and easy adjustment.
2. Using an Allen key (size 2 mm), increase the extension resistance by turning the valve screw (**Fig. 8, E**) clockwise until a comparatively long step is taken on the contralateral side and a harmonious gait pattern is achieved. For faster swinging forward on the prosthesis side, the valve screw is turned anticlockwise. Clicks and the scale on the valve screw enable precise and easy adjustment.

Never set the extension damping so high that full extension of the knee joint is no longer possible. The flexional joint position under load will inevitably cause the patient to fall!

The screws (**E**) and (**F**) must not be turned beyond the scale markings (**Fig. 8**).

■ Activating the high-resistance mode

To prevent unintended flexion of the joint (e.g. when moving through water), the joint has a high-resistance mode. This increases the flexion resistance so that flexion is only possible with very high force. The extension damping is not affected by this.

The high-resistance mode is activated by flipping the lever on the front of the damper (Fig. 8, H).

The high-resistance mode becomes ineffective when the freewheel mode (cycling function) is activated. The freewheel mode overrides all other hydraulic settings.

■ Activating the freewheel mode (cycling function)

The freewheel mode (cycling function) is activated and deactivated by pressing the push-button on the side of the hydraulic damper.

Do not touch the joint mechanism when flexing the knee joint – **risk of injury from pinching fingers!**

Adjust the walking mode using the push-button (Fig. 9, G).

Adjust the freewheel mode using the push-button (Fig. 9, L).

■ Factory settings

Setting of swing phase control on delivery:

- On delivery, the joint is in active walking mode.
- On delivery, the high-resistance mode is deactivated.

Flexion: the valve screw is 6 clicks away from the open position (Fig. 8, F)

Extension: the valve screw is 0 clicks away from the open position (Fig. 8, E)

■ Maintenance and service

Generally, all moveable parts must be inspected visually and for functionality every 6 months; all non-moveable structural parts every 12 months. If the product is also used for sporting activities, the inspection intervals must be reduced to a maximum of 3 months. If the product is also used in water, the maintenance intervals must be reduced to a maximum of 3 months. We recommend the special joint grease lubricant for lubricating the gliding metal surfaces of the knee joint.

In the event of malfunctions, the product must be sent in for maintenance.

Please note that repairs to the product may be carried out only by Streifeneder ortho.production GmbH.

■ Declaration of conformity

Streifeneder ortho.production GmbH declares as manufacturer with sole responsibility that the product complies with Regulation (EU) 2017/745 for medical devices (MDR).

■ Disposal

The product and its packaging must be disposed of in accordance with local or national regulations.

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Juhime tähelepanu sellele, et toote remonti tohib teostada ainult ettevõtte Streifeneder ortho.production GmbH.